

SURE FOOT

EQUINE STABILITY PADS

FACT SHEET - **HARD SLANT**

Holds the Angle

Genuine SURE FOOT Equine HARD SLANT Pads are color-coded with Pink DuraFlex and stamped with the official SURE FOOT logo on one side, and ivory on the other side. They are packaged in pairs along with a Warranty Card registration form and QR code for immediate access to the SURE FOOT® Equine website - for on-line product registration and education.



HARD SLANT are made of the same material as the SURE FOOT HARD Pads, but cut on an angle. The HARD SLANT may be used on either side; however, the pad is more stable used colored side up due to the angle it produces in relation to the ground. HARD SLANT will give to heat and pressure, but have no lateral give. Therefore they will, for the most part, maintain the angle under the weight of the horse. They can be used in any direction: heel-high, heel-low, or laterally (supination or pronation).

Horses can stand anywhere along the angle of the HARD SLANT Pads, from the thin leading edge to the full height (up to 3") at the back.

Horses that are particularly unsure or nervous

benefit from using them this way as the horse can stand on the leading edge, which is very thin. This allows the horse an opportunity to experience the feeling without the full height of the SLANT Pads or the added height of the rectangular HARD pads.

Pads are suitable for use on any surface with horses that are shod, bare foot and/or wearing boots. Pads will show wear over time. Nicks, tears and cuts are considered normal wear and tear and do not limit the effects of the pads in any way.

- **Benefit:** Long-term use to aid in improvement of overall horse health and well-being.
- **Color Coded:** Coated with Pink DuraFlex, a highly durable, anti-microbial coating.
- **Shape:** Rectangular shape cut on an angle along the length.
- **Size:** 10"x12"x3" (25cmx30cmx7.6cm) so that there is sufficient area for the horse's hoof divided into two angled pads.
- **Pad Weight:** Each pad weight 11.4 oz (326 grams).
- **Weight load limit:** Can be used with horses of any size, weight or breed.
- **Safety:** Pads are constructed of material that will not cause injury to the horse. However, failure to follow instructions may cause harm to both horse and human. Always follow the safety guidelines on the back of the insert sheet, website and/or DVD and instructional videos.
- **Instructions for use:** Available on www.surefootequine.com and on YouTube: SURE FOOT Equine Channel. DVD available upon request.
- **Breed or discipline:** May be used with any breed that has a hoof; horse, donkey or mule.
- **Serviceability:** While the pads are not impervious to damage (tears and cuts) they will continue to service the horse even when they look ugly due to use.
- **Portable:** Bring your pads to the arena, competition, field, veterinary clinic or anywhere you horse needs comfort and stability. The pads are easily used in any location.
- **General cleaning Instructions:** Easily cleaned with a soft brush, soap and warm water. Avoid scouring powder or other abrasive cleaners as they may damage the coated side.
- **Health Safety:** Disinfect with soap and water.

- **Hoof Preparation:** Bare foot horses tend to cut the pads so be sure to file off any sharp hoof edges to minimize damage. Shod horses tend to tear the edges of the pad if they slip off.
- **Uses:** Pads are beneficial to the horse pre-exercise, post-exercise, during a training session to improve and enhance performance, recovery from injury while stall bound or in limited movement, reset proprioceptors to establish normal patterns of movement, calming for any breed horse in any discipline, improve overall well-being of horse and response to instruction, improve relationship between horse and owner as the owner is recognized by the horse as caregiver.
- **Stackable.** SURE FOOT HARD SLANT Pads can be combined to stack layers generally two high but in some instances three high. Typically the HARD SLANT is placed as the top layer. Stacking allows the person to give the horse an experience that matches the different habit patterns of the body. Do not force a horse to stand on stacked pads. Stacked pads may slip as the horse steps off or slides the pads and this may cause additional wear and tear, therefore if you are stacking expect to see more damage. However, stacking allows the horse to experience limb differential and in many cases is beneficial to establishing a more normal, balance posture and movement.
- **Durability:** Genuine SURE FOOT pads are durable and have been designed and tested to hold up under the demands of use with horses in any environment. Pads from the original production run in January 2016 are still in use today.
- **Storage** – Be mindful of where you leave your pads. In very hot temperatures Hard pads will soften considerably. Do not leave any SURE FOOT pads in an enclosed car, horse trailer or other confined spaces in hot and/or humid temperatures. Doing so may void the warranty. A rule of thumb is not to leave your pads exposed to direct sunlight for prolonged a period of time. If you wouldn't leave your dog in that location, then don't leave your pads there either!
- **Warranty:** Warranted against delamination of the DuraFlex coating and/or pad layers. Excessive use with multiple horses on a daily basis will cause significantly more damage than the individual horse owner. However, the pads will continue to perform.
- **Manufacturing:** Made of as much environmentally friendly material as possible, designed and tested by Wendy Murdoch, carry the real SURE FOOT Logo and are assembled in the USA under the guidance and direction of The Murdoch Method, LLC.
- **Prop 65 Statement:** For information regarding the State of California Proposition 65 please see the description of the Pads at murdochmethod.com.
- **Warning:** *Please be sure to familiarize yourself with the SURE FOOT® Equine Stability Program before using SURE FOOT pads with your horse. Failure to follow instructions could result in death or serious injury. Check with your veterinarian before using SURE FOOT for any health-related conditions. Keep your hands away from the horse's hoof. Use your foot to position the pads. See instructions for use on the SURE FOOT DVD, SURE FOOT Equine website surefootequine.com and Sure Foot Equine [YouTube](#) Channel.*